



1 syringe



Eye Protocol

0.5cc per point


1+1 syringe



GENERAL Protocol

0.2cc per point

1+1+1 syringe



Enhancing Protocol

Inject 0.5cc per point around the eyes
Inject 0.25cc per point on the other points


HOW TO INJECT GOURI

General Protocol

☉ Inject to deep dermis

☉ Inject 0.2 cc per point


- **1st point:** Middle point of the forehead, above the center of an eye pupil
*On the forehead area, we recommend rubbing gently on the injection area after injection for 5-10 min
- **2nd point:** Horizontal with nostrils, about 2-3 cm away from the nose on the medial cheek
- **3rd point:** Same level as 2nd point, below the cheekbone
- **4th point:** Same level with the end of the lips line on the superior jowl, located in the middle between the 2nd and 3rd point
- **5th point:** Inject one point within the red marked areas, about 1.5-2.0 cm horizontally and vertically away from the of the eyes.



HOW TO INJECT GOURI

Eye Protocol

- Inject to deep dermis
- Inject 0.5 cc per point





- Injection point: Inject one point within the red marked area, about 1.5-2.0cm horizontally and vertically away from the of the eyes.

HOW TO INJECT GOURI

Enhancing Protocol

- ⦿ Inject to deep dermis
- ⦿ Inject 0.5 cc per point – eyes
- Inject 0.25 cc per point – entire face



- If you want special care, inject following the same points of General protocol with 0.25cc for 1st - 4th points and with 0.5cc for eyes